



**Minimize Imperfections,  
Naturally Improve**



**MINI**  
Facial System by Clareblend, Inc.

Device has a 1 year limited warranty.

**YOU MUST REGISTER YOUR DEVICE  
TO ACTIVATE THE WARRANTY.**

Activate warranty at:  
[www.clareblend.info](http://www.clareblend.info)

Christine Byer Esthetics  
Support@ChristineByer.com  
317-562-0222

Items included with MINI Facial System:

- 1 - MINI travel bag
- 1 - MINI device
- 1 - 2 oz conductive gel
- 1 - AC adapter
- 1 - manual
- 1 year limited warranty

The MINI is supplied with a wall plug adapter meaning there is no battery to replace and no loss of power due to low battery. Just plug and go without setting adjustments.

The USB connector on the cord can also be plugged into the following USB ports:

- Cell phone wall adapter
- External battery packs
- Computer USB port
- USB car charger

USB Compatible Power Sources



Support@ChristineByer.com

Please send email to get support with any repair requests. Christine Byer Esthetics will coordinate a repair request for you.

WARNINGS & PRECAUTIONS:

The MINI facial procedure is intended for facial stimulation and should not be used if the following indications are present:

- Pregnancy
- Seizures
- Cancer or tumors
- Cardiac pacemaker, or any electronic implanted defibrillator, simulators, or electronic devices as possible interference with these devices may occur
- Do not use if prone to phlebitis, thrombosis, or on varicose veins

Use MINI only on normal skin. Do not use on skin that is swollen, infected, or has cancerous lesions.

MINI is designed for cosmetic use only and for individuals in good health. As with any medical condition or after facial surgery of any kind, please consult your physician before use.

MINI is for adult use only and should be kept out of the reach of children.

Wait a minimum of 10 days before using MINI after the following procedures:

- Botox
- Injections and/or fillers
- Micro-needling (use after skin has fully healed)

To clean, wipe with a damp cloth.

DO NOT SUBMERGE IN WATER.

Troubleshooting MINI:

If no tone is heard when touching both electrodes to the skin:

1. Apply more gel
2. Ensure plug is secure in the adapter and power source
3. Ensure both electrode rods are touching the skin as both metal rods must be in contact with the skin for current to flow

When using the MINI, current should be sub-sensory, though some people may experience a slight tingling sensation. Protection circuitry protects the user and internal circuits from surges of alternating current or power source failure.

No adverse effects have been reported from MINI use. Mild redness from the stimulation of circulation may occur. If irritation or burning is experienced at any time, discontinue use of the device.

Clareblend conductive gel is part of the skin enhancement process. Read all gel ingredients prior to application. Do not use conductive gel if any ingredient(s) may trigger a known allergic reaction. For a patch test, apply a small amount of gel on the skin for 15 minutes or longer. If allergic reaction occurs or develops after 24 hours, discontinue use and contact your skin care professional for alternative product suggestions. Apply gel as needed to maintain smooth glide and conductivity.

Any harm, to any area of the body or face, resulting from the misuse or application, because of improper voltage settings, contamination of conductive gel, lack of cleaning electrodes, or other improper applications, are not the responsibility of Clareblend or any of its affiliates.

When procedure is complete, unplug device from power source, clean with damp cloth, and store in the provided travel bag. Apply a minimum SPF 15 moisturizer if going outside. Contact your skin care professional about supportive skin care products.

Device has a 1 year limited warranty. Activate warranty by submitting the online form found at:

<http://clareblend.info/mini-warranty-registration/> Warranty is void when device:

- Is transferred to another person
- Has been misused or subjected to accident
- Has improper maintenance or repair
- Has undergone an unauthorized modification

Cleanse the entire face and neck area. Apply a thin layer of conductive gel to one side of the neck & jaw line. Touch MINI electrodes to the skin at the base of the neck & wait for a single beep. Slowly glide upward toward jaw line. Maintain electrode contact with the skin while lifting & holding at the end of the movement until a double beep is heard. Move to the next section & continue movements, as shown, 3 times each.

Avoid thyroid and the center area of the neck as indicated as a rectangle on the diagram.



Movement 1

Position MINI near the center of the neck and touch the electrodes to the skin. Wait for a single beep, then move the device slowly toward the side of the neck. Lift & hold until a double beep is heard. Continue movements as shown, working each section 3 times on the neck area.

Avoid thyroid and the center area of the neck as indicated as a rectangle on the diagram.



Movement 2

Apply a thin layer of conductive gel to the jawline, cheek and temple area. Position MINI at the jaw line and touch both electrodes to the skin. Slowly glide the MINI up and under the cheek bone. Lift and hold until a double tone is heard. Move to the next section as shown by the arrows in the diagram and repeat process. Continue movements as shown, working each section 3 times.



Movement 3

Touch electrodes to the skin near the jaw line and slowly glide across the skin toward the side of the face. Bottom movement follows the jaw line with one electrode above & one electrode below the jaw line. Maintain contact with the skin to lift & hold at the end of the movement until a double beep is heard. Move to the next section as shown by the diagram and repeat until each movement has been completed 3 times.

It is okay to overlap movements.



Movement 4

Gently apply a thin layer of conductive gel around the eye, temple and brow area. Touch electrodes to the skin near the nose. After a single beep is heard, move MINI toward the ear and continue movement gliding upward treating crows feet area. Hold and lift at the end of the movement until a double beep is heard. Repeat 2 more times.

Use lighter pressure under the eye.



Movement 5

Touch electrodes to the skin at the base of the brow bone. Wait for a single beep, then lift at the brow and hold until a double beep is heard. Repeat three times in each section as arrows show on the diagram.

All eye movements should be on the orbital bone (never on the eyelid).



Movement 6



Apply a thin layer of conductive gel to one side of the forehead. Touch electrodes to the skin between the brows and wait for a single beep. Begin movement between the brows, moving up and over brow toward the side of the forehead. At the end of the movement, hold the electrodes in place against the skin while lifting until a double beep is heard. Movements for each section should be completed 3 times.



Movement 7

Touch electrodes to the skin at the top of the brow and wait for a single beep. Begin movement at the top of the brow and glide up toward the hairline. Lift & hold at the end of the movement until a double beep is heard. Complete each movement 3 times. After completing movements on one side of the face, move to the next side by applying a thin layer of conductive gel then begin movements 1-8. When finished with all movements on both sides of the face, apply sun protectant, unplug MINI from power source, clean electrodes with a damp cloth & store MINI in its travel bag.



Movement 8